



醫院管理局
HOSPITAL
AUTHORITY



Hospital Authority Convention 2014

Speed Presentation (SPP-P6.18)

Teddy LEE EO I, CSSD, TMH

7 May 2014





醫院管理局
HOSPITAL
AUTHORITY



Implementation of a Quality Management System in Sterile Supply Service is Critical to Continuous Quality and Safety Improvement

*Samuel Law TH, Lee CK, Wan CH, Fong YL, Li KY,
Yeung WC, Lam PI, Yip WY, Yip SY, Wong NP,
Wong WY*



*Central Sterile Supplies Department, New
Territories West Cluster*



 屯門醫院
TUEN MUN HOSPITAL
醫院管理局新界西醫院聯網

小欖醫院
SIU LAM HOSPITAL

Introduction



Quality Management Systems
for Medical Devices
ISO 13485:2003

Reference to
Policy of
Corporate

Quality Manual
Policies and procedures

Decontamination
Equipment

Management

- Washer Disinfectors
ISO 15883
HTM2030 (UK)
 - Sterilizers
EN 285, ISO 17665
HTM2010 (UK)
- Other Equipment
ISO 14937


The Facility

- Design
HBN13 (UK)
- Environment
ISO 14644

Risk Management
Training
Resources Monitoring
Auditing
Customer Focus
Tracking & Tracing of
Instruments

Quality Manual

1. Organization profile
2. Quality Policy and Objective
3. Management Responsibility
4. Human Resources Management
5. Material Management System
6. Product Standard
7. Risk Management, Measurement, analysis & Continue Quality Improvement
8. Document Control

	Cluster Central Sterile Supply Department	Document No.	NTWCSDD-A-QM-001-V5
	Quality Manual of Cluster Central Sterile Supplies Department	Version:	5
		Effective Date	01 APR 2014
		Next Review	01 APR 2016
		Page	Page 1 of 4

Contents

- 1.0 Introduction and Department Profile
 - 1.1 Introduction
 - 1.2 Department Profile
 - 1.21 Location
 - 1.22 Scope of Service
- 2.0 Quality Policy and Objective
- 3.0 Management Responsibility
 - 3.1 Management Commitment
 - 3.11 Mission
 - 3.2 Organization Chart
 - 3.3 Customer Focus
 - 3.31 External & internal communication
 - 3.32 Nearly Missed Event, complaint, feedback management
 - 3.33 Regular meeting
 - 3.4 Management audit and review
 - 3.41 Audit and review on management system
 - 3.42 Review on SOP
 - 3.43 Review on instructions & guidelines
 - 3.44 Review on customer information & feedback
 - 3.45 Improvement & Recommendations
- 4.0 Human Resources Management
 - 4.1 Job description
 - 4.11 administrative staff
 - 4.12 frontline staff
 - 4.2 Duty Arrangement
 - 4.21 Duty Hours
 - 4.211 administrative Staff
 - 4.212 frontline staff
 - 4.22 Duty Request Guideline
 - 4.23 Overtime guideline
 - 4.3 Leave Arrangement

Nearly Missed Case Management



	Cluster Central Sterile Supply Department	文件號碼: NTWCSSD-C-OM-003-V1
	Quality Improvement on Reporting Incident in NTWC CSSD	版本: 1
		生效日期: 01 APR 2014
		下次修訂: 01 APR 2015

Quality Improvement on Reporting Incident in NTWC CSSD



All incidents should be reported thro' email within 48 hours .

Severity	Incident belongs to group 3 & 4 should report Ms. Li and Mr. Law verbally.
1. very low risk	Incident stopped before reach to end-user.
2. Low risk	Incident related to non-conformity
3. Moderate risk	Incident related to expensive instruments / potential complaint
4. High risk	Incident may lead to delay in OT service/patient safety



Date	Incident Code	Finish Good Code	Washed by	Packer	Checker	Incident Details	Action Taken	Picture	Recorded by
03/03/2011	P11003/T1017			Dhan Mng	Tai Wong Ho	Foreign Body and dirt was found in a saw & air hose of Sterile Tracking Set T3017... At around 11:30 today clerk received call from the name of Theatre 8, who complained that there was dirt on the halls of an air hose of tracking set Hall's Reciprocating Saw. Then clerk immediately took another same type tracking set for replacement. Upon opening replacement set, it was also found that foreign bodies (whitish substance) on the jaw of the saw.	1. Root Cause- Packers are careless - 2. Immediate Action- Another same type tracking set for replacement - 3. Recommendation- All packers were informed of the incident and the procedures of checking washing the air hose and checking the saw after washing were reinforced - 4. Final Review (Date and Review)- The case would be included as training material for new comers - 3/6/2011 by EON Man Chan		Mary...

博愛醫院附屬中心消毒物品及器械物品供應部	文件號碼: OTWCSSD-C-OM-003-V1
包裝工作注意事項	版本: 1
	生效日期: 01 APR 2014
	下次修訂: 01 APR 2015

Incident Code: P12005D Finish Good Code: T1033
所有包裝員及審核員須徹底檢查手術儀器是否潔淨。凹位必須徹底清潔及吹水，留意有沒有骨粉




凹位必須徹底清潔及吹水



Quality Improvement Manual 品質提升手冊
Incident Code: T1003D Finish Good Code: T3017 所有儀器必須徹底檢查及清潔

Checklist of OT Set	Page: 1	
OT Set Code: T1004 (93)	Master Code: C 192	
	Washer ID: WD 4	
OT Set Name: ORTHO MINOR SET		
11009D: MINI CLIP 發現有異物，包裝時必須徹底檢查及清潔。		
12004D: 所有包裝員在使用時，必須注意份量，避免過多。		
P1203D: 所有包裝員在使用時，必須注意份量，避免過多。		
Section	Qty	Instrument Description
1	1	GRIP HOLDING, SPONGE, HANDLE, 250MM, 3/8"
(PRE-OPEN)	1	CLIP, MINI, 'D'
	1	CLIP, MINI, 'L'
	2	HANDLE, R. P., WD. 3, 88073R/PRO3/110-03/04, 103, 00/04
	1	-1000
	4	HANDLE, R. P., WD. 4, 88084R/802982/04, 104, 00/0P, 100-0
	1	W
2	4	FORCEPS, DISSECTING, ADDON, MICRO, 'T', 80078

Key Performance Indicator





醫院管理局
HOSPITAL
AUTHORITY



新界西·醫院聯網
New Territories West Cluster

Thank you!



 屯門醫院
TUEN MUN HOSPITAL
醫院管理局新界西醫院聯網

小欖醫院
SIU LAM HOSPITAL