#### Reduction of Allogeneic Blood Transfusion in Cardiac Surgery by Lowering Cardiopulmonary Bypass Prime Volume



Ng Wellon RN, Perfusion Team Division of Cardiothoracic Surgery Department of Surgery Prince of Wales Hospital

# Introduction



## **Cardiopulmonary Bypass (CPB)**



**Arterial Filter** 



Oxygenator

#### **Pump Head**







Panel



# Objective





- \* To minimise allogeneic blood transfusion for small size adult in priming and during CPB
- \* To maintain adequate organ perfusion during CPB, i.e. haematocrit (HCT) kept between 24-30%.

# Methodology



# Methodology

	Conventional group	Low prime group
Period	Mid-2015 to Mid-2016	Mid-2016 to November 2017
No. of patient	316	285
Body surface area (BSA)	< 1.8	
Arterial filter	Not integrated	Integrated
Priming volume	<b>1400</b> ml	1000 ml

### **Results**



### **Results**

No significant differences between both groups \*Baseline characteristics \*Body surface area (BSA) \*Type and urgency of the operation \*Perfusion techniques \*Haemotologic profiles

### **Blood Transfusion Requirement in Priming Cardiopulmonary Bypass Circuit**





#### Blood Transfusion Requirement during Cardiopulmonary Bypass





### Blood Transfusion Requirement after the Operation

# They were similar in both Conventional and Low Prime group

- \* First 24 hour Haematocrit
- **\* Discharge Haematocrit**

# Conclusion

![](_page_13_Picture_1.jpeg)

#### **Indications for Blood Transfusion during CPB**

![](_page_14_Figure_1.jpeg)

## **Optimizing CPB Practice**

	Capiox FX15R	Capiox FX 25R	Medtronic Affinity Fusion
Membrane Type	Microporous Polypropylene Hollow Fiber		
<b>Arterial Filter</b>	Integrated Arterial Filter		
Static Priming Volume	<b>135</b> ml	250ml	260ml
Blood Flow Range	<b>0.5-5.0</b> L/min	0.5-7.0 L/mi	n 1.0-7.0 L/min

## **Optimizing CPB Practice**

![](_page_16_Figure_1.jpeg)

**Small oxygenator** 

![](_page_16_Picture_3.jpeg)

**BSA > 1.8** 

Large oxygenator

### References

- Gravlee GP, Davis RF, Hammon H, & Kussman B (2016). Cardiopulmonary Bypass and Mechanical Support: Principles and Practice (4th ed.). Philadelphia, PA: Wolters Kluwer Health/Lippincott Williams & Wilkins
- Lahanas A, Argerakis PW, Johnson KA, Burdan ML & Ozdirik JE. A retrospective comparison of blood transfusion requirements during cardiopulmonary bypass with two different small adult oxygenators. Perfusion 2013; 28(6) 541-545.
- Ranucci M, Pistuddi V, Carboni G, et al. Effects of priming volume reduction on allogeneic red blood cell transfusions and renal outcome after heart surgery. Perfusion 2015; 30 (2) 120-126.
- Bronson SL, Riley BJ, Blessing JP, et al. Prescriptive patient extracorporeal circuit and oxygenator sizing reduces haemodilution and allogeneic blood product transfusion during adult cardiac surgery. The Journal of ExtraCorporeal Technology 2013; 45 167-172.

![](_page_18_Picture_0.jpeg)